

## 6 What to do if...



Error message,  
OFL

| Cause                                   | Remedy                             |
|---|------------------------------------|
| pH electrode:                           |                                    |
| – Not connected                         | – Connect electrode                |
| – Air bubbles in front of the diaphragm | – Remove air bubbles               |
| – Air in the diaphragm                  | – Extract air or moisten diaphragm |
| – Cable broken                          | – Replace electrode                |
| – Gel electrolyte dried out             | – Replace electrode                |

Error message,  
E3

| Cause  | Remedy                                       |
|--|--|
| pH electrode:  |  |
| – Diaphragm contaminated                                     | – Clean diaphragm                            |
| – Membrane contaminated                                      | – Clean membrane                             |
| – Moisture in the plug                                       | – Dry plug                                   |
| – Electrolyte obsolete                                       | – Replenish electrolyte or replace electrode |
| – Electrode obsolete   | – Replace electrode                          |
| – Electrode broken   | – Replace electrode                          |
| Measuring instrument:  |  |
| – Incorrect calibration procedure                            | – Select correct procedure                   |
| – Incorrect solution temperature (without temperature probe) | – Set up correct temperature                 |
| – Socket damp  | – Dry socket                                 |

|                                 |                              |   |
|---------------------------------|------------------------------|---|
|                                 | Buffer solutions:            |   |
|                                 | – Incorrect buffer solutions | – Change calibration procedure                    |
|                                 | – Buffer solutions too old   | – Only use once. Note the shelf life              |
|                                 | – Buffer solutions depleted  | – Change solutions                                |
| <b>No stable measured value</b> | <b>Cause</b>                 | <b>Remedy</b>                                     |
|                                 | pH electrode:                |   |
|                                 | – Diaphragm contaminated     | – Clean diaphragm                                 |
|                                 | – Membrane contaminated      | – Clean membrane                                  |
|                                 | Sample:                      |   |
|                                 | – pH value not stable        | – Measure with air excluded if necessary          |
|                                 | – Temperature not stable     | – Adjust temperature if necessary                 |
|                                 | Electrode + sample:          |   |
|                                 | – Conductivity too low       | – Use suitable electrode                          |
|                                 | – Temperature too high       | – Use suitable electrode                          |
| – Organic liquids               | – Use suitable electrode     |   |
| <b>LoBat</b>                    | <b>Cause</b>                 | <b>Remedy</b>                                     |
|                                 | – Batteries almost depleted  | – Replace batteries (see section 5.1 MAINTENANCE) |

|   |  |  |
|---|--|--|
| <b>Obviously incorrect measured values</b>                                | <b>Cause</b>   | <b>Remedy</b>  |
|   | pH electrode:  |  |
|   | – pH electrode unsuitable                                    | – Use suitable electrode   |
|   | – Temperature difference between buffer and sample too large | – Adjust temperature of buffers or samples   |
|   | – Measuring procedure not suitable                           | – Follow special procedure   |
| <b>Instrument does not react to keystroke</b>                             | <b>Cause</b>   | <b>Remedy</b>  |
|   | – Operating state undefined or EMC electric stress unallowed | – Processor reset:<br>Press the  key and switch on instrument                   |
| <b>You would like to know which software version is in the instrument</b> | <b>Cause</b>   | <b>Remedy</b>  |
|   | – e.g. question of the WTW service department                | – Press the  key and switch on instrument. The software version is displayed. |