
Operating manual
Bedienungsanleitung
Mode d'emploi
Instrucciones de operación

pH 330i pH 340i

pH/mV Pocket Meter	Page	3
pH/mV Taschenmeßgerät	Seite	17
pH/ mV- mètres de poche	Page	31
pH/ mV Medidor manual	Página	45

**Note to this
operating manual**

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pH 330i/340i - List of contents

Safety	4
Display and jack field	4
Mains power (pH 340i only)	5
Switching on the measuring instrument	5
Measuring	6
Calibration	7
Saving data	11
Outputting the data memory	11
Transmitting data (only pH 340i)	11
Configuring	12
Resetting to default settings	13
Technical data	14

Safety

Safety instructions


The individual chapter of this operating manual use the following safety instruction to indicate various types of danger:

Warning


indicates instructions that must be followed precisely in order to avoid the possibility of slight injuries or damage to the instrument or the environment.

Display and jack field

Display

Jack field

- | | |
|--|--|
| 1 pH electrode | 4 Temperature sensor |
| 2 Reference electrode | 5 RS232 interface or analog output/
recorder (only pH 340i) |
| 3 Plug-in power supply (only
pH 340i) | |

Warning

Only connect electrodes to the measuring instrument that do not return any unallowed voltages or currents (> SELV and > current circuit with current limiting).

Almost all electrodes - in particular WTW electrodes - fulfill these conditions.

Mains power (pH 340i only)


You can either operate the measuring instrument with batteries or with the plug-in power supply. The plug-in power supply provides the measuring instrument with low voltage (7.5 V ... 12 V DC). This saves the batteries.

Warning

The line voltage at the operating site must lie within the input voltage range of the original plug-in power supply (see TECHNICAL DATA).

Warning

Only use original plug-in power supplies (see TECHNICAL DATA).

- Plug the jack (1) into the socket (2) of the measuring instrument.
- Connect the original WTW plug-in power supply (3) to an easily accessible mains socket.

Switching on the measuring instrument

- Press the <  > key.
Display test appears briefly on the display.
After this, the measuring instrument automatically switches to the measuring mode. The display shows the relevant measured value.

Measuring

Overview of the measuring modes:

Special functions:

Select pH resolution

If you press **<M>** while pressing the **<RUN/ENTER>** key, the resolution swaps between 0.01 and 0.001.

AutoRead (drift control)

The AutoRead function checks the stability of the measurement signal (except for the measurement of the ORP voltage). Activate AutoRead with **<AR>**. Press **<RUN/ENTER>** to start the AutoRead measurement. During the AutoRead measurement, *AR* flashes on the display until a stable measured value is reached. This can be terminated at any time taking over the current measured value with **<RUN/ENTER>**.

Temperature measurement during pH measurements

You can use pH electrodes with and without temperature sensors. The measuring instrument recognizes which temperature sensors are connected and switches automatically to the correct mode for the temperature measurement (display *TP*).

Manual temperature input: Set the temperature value with **<▲>** **<▼>**.

Calibration

AutoCal TEC

This process is specially adapted to the WTW technical buffer solutions as a fully automatic **single, two-point or three-point calibration**. The buffer solutions are automatically recognized by the measuring instrument. Depending on the instrument setting, the instrument displays the relevant buffer nominal value or the current electrode voltage in mV.

Valid buffer (values at 25 °C): 2.00 / 4.01 / 7.00 / 10.01

Note

The calibration for pH 10.01 is optimized for the WTW technical buffer solution TEP 10 Trace or TPL 10 Trace. Other buffer solutions can lead to an erroneous calibration. The correct buffer solutions are given in the WTW catalog or in the Internet.

AutoCal DIN

as AutoCal TEC, but matched to buffer solutions according to DIN 19266.

Valid buffer (values at 25 °C): 1.679 / 4.006 / 6.865 / 9.180

ConCal

is the conventional **two-point calibration** using two buffer solutions (pH 7.0 ± 0.5 and any other buffer solution) or a **single-point calibration** using any buffer solution and is used as a high-speed method.

Calibration evaluation

After calibrating, the measuring instrument automatically evaluates the calibration. The asymmetry and slope are evaluated separately. The worst evaluation appears on the display.

Display	Asymmetry[mV]	Slope [mV/pH]
	-15 ... +15	-60.5 ... -58
	-20 ... +20	-58 ... -57
	-25 ... +25	-61 ... -60.5 or -57 ... -56
	-30 ... +30	-62 ... -61 or -56 ... -50
Clean the electrode according to the electrode operating manual		
 Invalid calibration	< -30 or > 30	< -62 or > -50

**Note**

If a printer is connected to the interface, a calibration protocol is automatically printed after a valid calibration.

AutoCal TEC / AutoCal DIN procedure:

The following example shows the AutoCal TEC calibration. In the AutoCal DIN calibration, *Cd...* instead of *Ct...* appears on the display. Otherwise, both procedures are identical.


- Connect the pH electrode to the measuring instrument.
- Press the **<CAL>** key repeatedly until the display *Ct1* and the function display *AutoCal TEC* appears. The sensor symbol displays the evaluation of the last calibration (or no sensor symbol in the delivery state or after the measurement parameter has been reset).
- Immerse the pH electrode into the first buffer solution.
- When measuring without a temperature sensor:
Enter the temperature of the first buffer solution with **<▲>** **<▼>**.
- Press **<RUN/ENTER>**. The AutoRead measurement begins. The electrode voltage (mV) or the buffer nominal value appears on the display. If the measured value is stable, *Ct2* appears.

**Note**


At this point, the AutoCal TEC calibration can be terminated with **<M>**. This corresponds to a **single-point calibration**. To do this, the instrument uses the Nernst slope (-59.16 mV/pH at 25 °C) and determines the asymmetry of the electrode.

- Immerse the **two-point calibration** pH electrode in the second buffer solution to continue.
- When measuring without a temperature sensor:
Enter the temperature of the second buffer solution with **<▲>** **<▼>**.
- Press **<RUN/ENTER>**. The AutoRead measurement begins. The electrode voltage (mV) or the buffer nominal value appears on the display.
- If the measured value is stable, the instrument displays the value of the slope and the evaluation of the two-point calibration.
- Press **<RUN/ENTER>**. The instrument displays the value of the asymmetry.

**Note**

At this point, the AutoCal TEC calibration can be terminated with **<M>**. This corresponds to a **two-point calibration**.

- Immerse the pH electrode in the second buffer solution to continue the **two-point calibration**.
- When measuring without a temperature sensor:
Enter the temperature of the third buffer solution with **<▲>** **<▼>**.
- Press **<RUN/ENTER>**. The AutoRead measurement begins. The electrode voltage (mV) or the buffer nominal value appears on the display.
- If the measured value is stable, the instrument displays the value of the slope and the evaluation of the three-point calibration.
- Press **<RUN/ENTER>**. The instrument displays the value of the asymmetry.
- Switch to the measuring mode with **<M>**.

ConCal procedure:

- Connect the pH electrode to the measuring instrument.
- Press the <CAL> key repeatedly until the display *ASY* and the function display *ConCal* appears. The sensor symbol displays the evaluation of the last calibration (or no sensor symbol in the delivery state or after the measurement parameter has been reset).
- Immerse the pH electrode in the first buffer solution (pH 7.0 ± 0.5 in two-point calibration).
- When measuring without a temperature sensor:
Enter the temperature of the first buffer solution with <▲> <▼> while pressing the <RUN/ENTER> key.
- Press the <RUN/ENTER> key.
The measured pH value appears on the display.
- Set the measured value to the nominal pH value of the buffer solution (at the current temperature) with <▲> <▼>.
- If the measured value is stable, press the <RUN/ENTER> key. The value of the asymmetry appears.
- Press the <RUN/ENTER> key. *SLO* appears.

**Note**

At this point, the ConCal calibration can be terminated with <M>. This corresponds to a **single-point calibration**. To do this, the instrument uses the Nernst slope (-59.16 mV/pH at 25 °C) and determines the asymmetry of the electrode.

- Immerse the **two-point calibration** pH electrode in the second buffer solution to continue.
- When measuring without a temperature sensor:
Enter the temperature of the second buffer solution with <▲> <▼> while pressing the <RUN/ENTER> key.
- Press <RUN/ENTER>.
- Set the measured value to the nominal pH value of the buffer solution (at the current temperature) with <▲> <▼>.
- When the measured value is stable, press the <RUN/ENTER> key. The instrument displays the value of the slope and the calibration evaluation.
- Press <RUN/ENTER>. The instrument displays the value of the asymmetry.
- Switch to the measuring mode with <M>.

Saving data

Saving data manually

Press the **<STO>** key in the measuring mode (display *No.* with the number of the next free memory location). Then, press **<RUN/ENTER>** and enter the ID number with **<▲>** **<▼>**. Terminate the save with **<RUN/ENTER>**. The instrument changes to the measuring mode.

Saving data automatically *Int 1*

The portable measuring instrument pH 330i/340i has an internal data memory. It can store up to 800 datasets. The save interval (*Int 1*) determines the chronological interval between automatic save processes.

To set up the save interval, press **<STO>** while pressing the **<RUN/ENTER>** key (display *Int 1*) and set the interval with **<▲>** **<▼>** key. Then, press **<RUN/ENTER>** and enter the ID number with **<▲>** **<▼>**. Change to the last active measuring mode with **<RUN/ENTER>**. The automatic save is switched on (display *Auto Store*).

Outputting the data memory

You can output the data memory with the **<RCL>** key. By repeatedly pressing the **<RCL>** key, you reach the following functions:

- StO dISP* Output stored data on the display
- StO SER* Output stored data via the serial interface (only pH 340i)
- CAL dISP* Output calibration data on the display
- CAL SER* Output calibration data via the serial interface (only pH 340i)

Start the output with **<RUN/ENTER>**.

Transmitting data (only pH 340i)

Transmitting data manually

Press the **<RUN/ENTER>** key in the measuring mode. This manually triggers a data transmission of the current measured value to the serial interface at any time - independently of the selected intervals.

Transmitting data automatically *Int 2*

The interval to the data transmission (*Int 2*) determines the chronological interval between automatic data transmissions. After the selected interval expires, the current data record is transmitted to the interface. To set up the transmission interval, press **<RCL>** while holding down the **<RUN/ENTER>** key (display *Int 2*). Then, set the interval with **<▲>** **<▼>**.

Configuring

Note

You can leave the configuration menu at any time with **<M>**. The parameters already changed are stored.

Procedure for configuring (factory settings appear in bold typeface):

- Switch off the measuring instrument.
- While pressing the **<M>** key, press **<O>**.
- Select the required Baud rate with **<▲>** **<▼>**.
Selection: 1200, 2400, **4800**, 9600 Baud.
- Then press the **<RUN/ENTER>** key.
- Select the required display during the pH calibration with **<▲>** **<▼>**.
Selection: **Buffer nominal value (/pH)** or current electrode voltage (*mV*).
- Then press the **<RUN/ENTER>** key.
- Set the required calibration interval with **<▲>** **<▼>**.
Selection: 1... 7... **999** d.
- Then press the **<RUN/ENTER>** key.
- Select the date and time step-by-step with **<▲>** **<▼>**. In between, press the **<RUN/ENTER>** key each time.
- When the last **<RUN/ENTER>** key has been pressed, the instrument changes to the last active measuring mode.

Resetting to default settings

You can reset the measurement parameters and the configuration to the delivery status separately from one another (initialization).

Measurement parameters

The following measurement parameters can be reset to the delivery status:

pH/ORP (pH *Int*)


Measuring mode	pH
Slope	-59.16 mV/pH
Asymmetry	0 mV
Manual temperature input	25 °C
Resolution (pH display)	0.01

Configuration parameters

The following configuration parameters (*Int*) can be reset to the delivery status:

Baud rate	4800
Display during pH calibration	Buffer nominal value
Interval 1 (automatically saved)	OFF
Interval 2 (for data transmission)	OFF

Procedure of the reset:

- To switch it off, press <CAL> while pressing the <RUN/ENTER> key.
- To reset the measurement parameters with <▲> <▼>, select YES and confirm with <RUN/ENTER>, or
- Continue to the configuration parameters without resetting (*Int*) with the <RUN/ENTER> key.
- After the configuration parameters, *Int* changes the instrument to the last active measurement mode.

Technical data

Dimensions and weight	Length [mm]	172	
	Width [mm]	80	
	Height [mm]	37	
	Weight [kg]	approx. 0.3	
Mechanical structure	Type of protection	IP 66	
Electrical safety	Protective class	III	
Test certificates	cETLus, CE		
Ambient conditions	Storage	- 25 °C ... + 65 °C	
	Operation	-10 °C ... + 55 °C	
	Climatic class	2	
pH/ORP measuring ranges		Measuring range	Resolution
	pH	- 2.000 ... + 19.999	0.001
		- 2.00 ... + 19.99	0.01
	U [mV]	- 999.9 ... + 999.9	0.1
- 1999 ... + 1999		1	
T [°C]	- 5.0 ... + 105.0	0.1	
Precision of pH/ORP (± 1 digit)	pH (± 2 pH units from the calibration point)	± 0.003 at + 15 °C ... + 35 °C	
		± 0.01	
	U [mV]	± 0.2 at + 15 °C ... + 35 °C	
T [°C]	NTC 30: ± 0.1		
	PT 1000: ± 0.5 at 0 °C ... 15 °C ± 0.1 at 15 °C ... 35 °C ± 1 at 35 °C ... 55 °C		

pH correction function	Temperature input	
	Manually [°C]	- 20 ... + 130
Analog output (pH 340i only)	Automatic switchover when the recorder is connected by the cable, AK 323.	
	Output signal	corresponds to the electrode voltage, i.e. 50-62 mV/pH unit
	Accuracy	± 0.1 % of the measured value ± 4 mV
	Internal resistance	< 5 Ohm (current limited to max. 0.2 mA output current)
Serial interface (only pH 340i)	Automatic switchover when a PC or a printer is connected via the cable, AK 340/B or AK 325/S.	
	Type	RS232, data output
	Baud rate	Can be set to 1200, 2400, 4800, 9600 Baud
	Data bits	8
	Stop bit	2
	Parity	None
	Handshake	RTS/CTS + Xon/Xoff
	Cable length	Max. 15m

Power supply	Batteries	4 x 1.5 V alkali-manganese batteries, Type AA
	Operational life	approx. 3000 operating hours
	Power supply (only pH 340I)	The following applies to all plug-in power supplies: Connection max. Overvoltage category II Plug-in power supply unit (Euro, US, UK, Australian plug) FRIVO FW7555M/09, 15.1432 Friwo Part. No. 1822089 Input: 100 ... 240 V ~ / 50 ... 60 Hz / 400 mA Output: 9 V = / 1,5 A
Guidelines and norms used	EMC	EG guideline 89/336/EWG EN 61326 A1:1998 EN 61000-3-2 A14:2000 EN 61000-3-3:1995 FCC Class A
	Instrument safety	EG guideline 73/23/EWG EN 61010-1 A2:1995
	Climatic class	VDI/VDE 3540
	Type of protection	EN 60529:1991

FCC Class A Equipment Statement

Note: This equipment has been tested and found to comply with the limits for a Class A digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference when the equipment is operated in a commercial environment. This equipment generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instruction manual, may cause harmful interference to radio communications. Operation of this equipment in a residential area is likely to cause harmful interference in which case the user will be required to correct the interference at his own expense.

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Wissenschaftlich-Technische Werkstätten GmbH